

Sand Run Pharmacy

40 Sand Run Road
Akron, Ohio 44313

Phone: 330 - 864 - 2138 Website: www.sandunpharmacy.com

Topi-CLICK Micro

Veterinary PLO

User Guide

Cap



How to Use:

Step 1: Unscrew the cap at the top of the Topi-CLICK Micro® device.

Step 2: While holding upright (applicator up) turn the base until the device “*clicks*”. This noise indicates that a single dose has been dispensed.

- Each dose or *click* dispenses 0.05mL of the medication
 - 1 click = 0.05mL dose
 - 2 clicks = 0.1mL dose
 - 3 clicks = 0.15mL dose

Step 3: Apply to the hairless section of the ear. Continue rubbing **with the attached tip** until all medication is absorbed into the skin.

Step 4: Screw the cap back on to the top of the Topi-CLICK Micro® device.

*** Store the Topi-CLICK Micro® device at room temperature, do not expose to cold temperatures.



The Topi-CLICK Micro Device is recyclable

A safer way to deliver medication

The Topi-CLICK Micro® device allows the client to apply topical medications directly to the surface of the skin, *without* being exposed themselves.

About the Device:

- This product is intended for dispensing **topical medication ONLY**; do not administer any products orally from this device.
- The Topi-CLICK Micro® allows for multiple, accurately measured doses to be dispensed and applied with one device. The easy to use device replaces the practice of dividing medications into multiple syringes. The device is simple to use and delivers the same set dose with every use.
- The **REFILL Indicator** is marked on the device in order to ensure timely reordering of medication. When the black plunger reaches the line below the word “**refill**”, there are about 32 *clicks* left.



If you have any further questions/concerns while administering this product, contact our staff at Sand Run Pharmacy (330)-864-2138

~or~

Check out our videos on the Facebook Page “**Sand Run Pharmacy Prescription & Veterinary Compounding**” or on our website; www.sandrunpharmacy.com